

LASER THERAPY FOR HORSES

HORSEMIX OY

HAMINA
SUOMI, FINLAND

EFFECTS OF LASER LIGHT

- ANTI-INFLAMMATORY
- INCREASING RHEOLOGY OF BLOOD
- REDUCING EDEMA
- ANALGETIC
- STIMULATING OF BLOOD VESSELS PROLIFERATION



LASER THERAPY WAS USED IN HUMAN MEDICINE SINCE 1975

QUALITIES OF LASER LIGHT:

- PAINLESS
- ASEPTIC (NO CONTAMINATION)
- DEEP PENETRATING (5- 7 sm)
- WELL CONTROLLED (EASY TO COUNT THE DOSAGE)
- LOCALIZED (CAN BE APPLIED LOCALLY)

SEMICONDUCTOR LASER – HART OF THE SYSTEM


- COMPACT
- POWERFULL
- FLEXIBLE (CAN WORK WITH DIFFERENT TIPS)
- AUTONOMIC (CAN OPERATE IN STABLES)
- WATERPROOF



LASER SPECIALLY DESIGNED TO WORK IN STABLES: VERY SIMPLE CONTROLS



THERAPY LIGHT CHOSEN TO PENETRATE UP TO 10 sm DEEP



Picture 1.

ORIGINAL FIBER OPTICS FOR INNER ORGANS ACCSES



ALL FIBER OPTICS ARE FLRXIBLE AND PROTECTED FROM DUST AND MUD



FAVORITE 3000 IS FIXED IN ITS BAG ON THE THERAPIST'S CORPS



WE OFFER A CHOICE OF LASER HEADS TO FAVORITE 3000:

- 30 mW red laser for acupuncture



100 mW infrared laser module




Laser “SHOWER”, i.e cluster from 19 lasers joined in one module

- Infrared shower – 190 mW total power
- Red laser shower – 95 mW total power – better for surface treatments, such as skin sores, etc.



Laser shower, or DOUCHE can be equipped by fiber optic “comb”

- It serves to easier penetration of energy to the skin of animal
- The “comb” is interchangeable and easy to wash




Hi- power laser module 300 mW

- Infrared range for deep action
- Designed for treatment of large skin sores
- Ideal for treating joints and tendon injuries



Treatment is very easy and supported by original brochures



FAVORITE -3000 has a service support all over Finland



We are wishing You to get quick results of therapy by Your Favorite laser!

